

# ZOLL AED PLUS MAINTENANCE

**YEAR:**

Site Name:	Primary Liaison:
Site Address:	Backup Liaison:
Site Phone:	

SERIAL # OF AED:	Last time battery was changed (5 year battery life):
Adult Pad Expiration Date:	Pediatric Pad Expiration Date (optional)
Specific Location of Unit:	

Inspect unit at least once monthly and check for following.

- Verify the green check showing the unit is ready to use.
- Turn on AED to ensure it starts
- Verify electrodes are within their expiration date.
- Verify batteries are within their expiration date.
- Verify electrodes are pre-connected to the input connector.
- Verify supplies are available for use (razor, mask, gloves)

	January	February	March	April	May	June	July	August	September	October	November	December
<b>Check done</b>												
<b>Initials</b>												

**Comments:**

**KEEP THIS FORM FOR YOUR RECORDS**

**If any problems are found during inspection contact Zoll Technical Support at 1-800-348-9011  
For local support from Sedona Fire District contact Josh Wells, [padinfo@sedonfire.org](mailto:padinfo@sedonfire.org), 928-699-8853**